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

# Instructions for XXL-Rehab Minimaxx – Push motor

## CAUTION

This product should not be used without proper instruction from a healthcare professional.

Using the product for individuals weighing more than the maximum weight-bearing capacity may pose a risk for damaging the product.

Cobi Rehab assumes no responsibility for any damage or injury caused by improper installation, assembly or use of this product.

If components are damaged or missing, contact Cobi Rehab immediately. Do NOT use substitute parts.

Cobi Rehab cannot be held liable for faults or accidents occurring after repairs by service staff without written authorisation from Cobi Rehab.

 To ensure proper use of the XXL-Rehab Minimaxx Push motor, please read these instructions thoroughly before using the product. In doing so, pay particular attention to the warning notes marked with .

Also remember that in order to prevent the risk of injury to the user, the wheelchair may be operated only by trained caregivers or aides or by individuals who have been instructed in its proper use.

## USE

XXL-Rehab Minimaxx Push Motor is a PUSH motor which is designed for XXL-Rehab Minimaxx manual wheelchair. On the following pages you will find the user manual for the PUSH motor, referred to as PUSH motor.

### Component overview

On page 11 there is an overview of the components that have been delivered with PUSH motor.

### Assembly instructions

On page 12-15 there is an assembly instruction as to how PUSH motor is mounted on XXL-Rehab Minimaxx.

### Driving instructions

On page 16 there is an instruction as how to drive PUSH motor.

### Charging instructions

On page 17 there is an instruction as how to charge the battery of PUSH motor.

# Component overview

Showing all parts that have been delivered with PUSH motor

## Remote control


1. Remote control
2. Remote console
3. Thumbscrew
4. 2 pcs hardware bolts
5. Upper handle bracket
6. 3 pcs spacers
7. Speeder
8. Lower handle bracket

## Wheelchair bracket

9. 2 pcs hub bushing
10. 2 pcs locknuts for bush
11. 2 pcs splice
12. 2 pcs thumbscrew
13. 2 pcs protection caps

## Push motor

14. Cord for remote control
15. Battery console
16. Lifting strap
17. Pole
18. Profile tires
19. Battery
20. Battery charger

# Assembly instructions

Assembly instructions will demonstrate how PUSH motor is mounted correctly on XXL-Rehab Minimaxx.


Before use of PUSH motor it is important that the battery is charged.

See **CHARGING INSTRUCTIONS** on page 17.

**Fig. 1**

Start by mounting the battery (19) on the PUSH motor battery console (15). When the battery is pushed to the bottom of the console it is locked.


**Fig. 1**

**Fig. 2**

Before PUSH motor is mounted on the XXL-Rehab Minimaxx frame we recommend to lock the parking brakes on XXL-Rehab Minimaxx.


**Fig. 2**

**Fig. 3**

Mount the two hub bushes (9) and the locknuts (10) on XXL-Rehab Minimaxx. Screw the hub bushes in the front threaded hole on XXL-Rehab Minimaxx right and left frame. Mount the locknuts on the inside of the frames.

**Fig. 3**


# Assembly instructions

**Fig. 4**

First remove the protection caps (13). Then slide the splices (11) with thumbscrew (12) onto the pole (17). PUSH motor is placed underneath XXL-Rehab Minimaxx so the ends of the splices are in front of the hub bushes.


**Fig. 4**

**Fig. 5**

Slide the splices out to each side until the ends rest in the hub bushes. Lock the splices onto the pole with the thumbscrews. PUSH motor is now mounted on the XXL-Rehab Minimaxx frames.


**Fig. 5**

**Fig. 6**

Mount the PUSH motor remote control on XXL-Rehab Minimaxx right push handle (23). Push a suitable spacer into the push handle tailpipe. Upper and lower handle bracket is tightened around the spacer with the two supplied mounting bolts (4).

**Fig. 6**

# Assembly instructions

**Fig. 7**

Screw PUSH motor speeder on the lower handle bracket with the thumbscrew (3). Push the remote console (2) over the upper handle bracket. Click the remote control (1) on the console. Connect the connectors on the accelerator and remote control.


**Fig. 8**

The remote control is now assembled. With the right hand on the push handle it is possible to reach the speeder and at the same time have a good grip on XXL-Rehab Minimaxx.

**Fig. 9**

Connect PUSH motor remote control to PUSH motor. Stretch the spiral cord and connect it to the PUSH motor cord of the remote control. Remote control and PUSH motor is now ready for use.

**Fig. 7**

**Fig. 8**

**Fig. 9**

# Assembly instructions

**Fig. 10**

PUSH motor can be lifted from the driving surface when not needed. A lifting strap hook (21) is mounted on XXL-Rehab Minimaxx left push handle (24).

**Fig. 11**

The lifting strap (16) is mounted on PUSH motor. Adjust the length of the strap so the wheels on PUSH motor are lifted approx. 2-3 cm over the driving surface when the lifting strap hangs on the lifting strap hook.

**Fig. 12**

Release XXL-Rehab Minimaxx parking brakes and take the lifting strap of the lifting strap hook. Now PUSH motor is ready for use.

**For the best possible start we recommend to read DRIVING INSTRUCTIONS on page 16.**

**Fig. 10**

**Fig. 11**

**Fig. 12**

# Driving instructions

We recommend that the first ride with PUSH motor takes place in a location with plenty of room, smooth running surface and no passenger in the wheelchair.

**Fig. 1**

Check that the PUSH motor battery is fully charged before removing from the battery charger. If the charge lamp on the charger turns green, the battery is fully charged. However, if the charge lamp is red, charging is still underway and the drive may need to be considered.

**See CHARGING INSTRUCTIONS on page 17.**

Lower the PUSH motor onto the driving surface and position the lifting strap so it cannot be entangled in the wheels.

**Fig. 2**

Turn on the power on the PUSH motor remote control by using the ON/OFF button (1). The battery indicator of the remote control (2) shows with a dash of green light how much power is left on the battery. With the buttons (3) and (4) you determine whether PUSH motor should drive forward or backward. The blue lamps (5) and (6) indicate the selected direction.

Adjust the PUSH motor top speed on the speed regulator of the remote control (7). Select "turtle" for low speed and "hare" for high speed. We recommend starting with "turtle".


Carefully lift up the remote control speeder and PUSH motor begins to push the wheelchair. When you let go of the speeder the wheelchair stops instantly. Wanting more speed push the speed regulator to "hare". Going downhill let slightly go of the speed regulator to reduce speed. On very steep slopes the wheelchair must be driven backwards down the hill.

**⚠ WARNING**

**To avoid accidents please remember:**

- To turn off the ON/OFF button during longer run breaks. This avoids errors if the wheelchair is vacated.
- To lock the parking brakes of the wheelchair when the wheelchair is vacated.

**Fig. 1**

**Fig. 2**

# Charging instructions

We recommend that PUSH motor is being charged after each drive so the battery is fully charged and ready for the next drive.

**Fig. 1**

The battery charger fits into a standard 230V wall outlet and the charging connector on the battery.

**Fig. 2**


The light on the charge lamp is red while charging is in progress. When the battery is fully charged the charge lamp is green. Be aware that the charge lamp is also green if the charger gets 230V but is taken out of the battery.

With normal use of PUSH motor charging the battery will take 6-8 hours.  
If the battery is run almost empty, charging will take up to 12 hours.


**⚠ WARNING**

**To avoid accidents and power failures, please remember:**


- PUSH motor can drive while being charged.
- PUSH motor must be fully charged both before and after a longer period of not being used.
- PUSH motor use power when the remote control is ON.

**Fig. 1**

**Fig. 2**

**MANUFACTURER AND DEALER**

Cobi Rehab  
Fuglebækvej 1D  
DK-2770 Kastrup  
Denmark  
Tel.: +45 70 25 25 22  
www.cobirehab.com

04/2018/JP

# COBI REHAB<sup>®</sup>

The Positioning Company

**COBI REHAB**  
Fuglebækvej 1D  
DK-2770 Kastrup  
Denmark  
Tel. +45 7025 2522  
[www.cobirehab.com](http://www.cobirehab.com)